

UCSD Center for Epigenomics ATAC-seq protocol for cells (adapted from Buenrostro et al., 2015; PMID: 25559105)

1. Thaw cells at 37°C for 2 min. Transfer to a 1.5 ml tube and spin down for 5 min at 500 g at 4 °C.
2. Immediately resuspend in 250 µl nuclei permeabilization buffer (NPB) and rotate for 5 min at 4°C.

Composition of nuclei permeabilization buffer (NPB):

Reagent	Stock concentration	Final concentration	Amount per 10 ml
BSA (A7906, Sigma)	--	5 %	500 mg
IGEPAL-CA630 (I8896, Sigma)	10 % (m/V)	0.2 % (m/V)	200 µl
DTT (D9779, Sigma)	200 mM	1 mM	50 µl
cOmplete EDTA-free protease inhibitor (05056489001, Roche)	25x	1x	400 µl
PBS (10010-23, Thermo Fisher Scientific)	--	--	9.35 ml

3. Spin down for 5 min at 500 g at 4 °C
4. Resuspend nuclei in cold 25 µl tagmentation buffer (TB) and place on ice. Take an aliquot and determine nuclei concentration using a hemocytometer, e.g. use 5 µl nuclei suspension, add 20 µl TB and 25 µl Trypan Blue Stain (0.4 %, 15250061, Thermo Fisher Scientific)

Composition of tagmentation buffer (TB)

Reagent	Stock concentration	Final concentration	Volume per 10 ml
Tris-acetate (pH=7.8, BP-152, Thermo Fisher Scientific)	1 M	33 mM	330 µl
K-acetate (P5708, Sigma)	3 M	66 mM	220 µl
Mg-acetate (M2545, Sigma)	300 mM	11 mM	367 µl
DMF (DX1730, EMD Millipore)	--	16 %	1.6 ml
Molecular biology water (46000-CM, Corning)	--	--	7.483 ml

5. Adjust nuclei concentration to 2,500 nuclei/µl (50K total) and transfer 20 µl for tagmentation to a 1.5 ml tube.

Tagmentation and library generation

6. Add 1 μl Tn5 mix 5 times by pipetting and incubate 60 min with 500 rpm at 37 °C
7. After completion add 100 μl PB buffer (Qiagen) and 5 μl Na-acetate (3 M, pH = 5.2). Purify using MinElute PCR Purification Kit (28004, Qiagen) and elute in 10 μl EB.
8. Amplify fragments by PCR. Please see "Appendix" for primer sequences.

Reaction setup

Reagent	Volume per sample
Tagmented DNA	10 μl
NEBNext 2x PCR MasterMix (M0541, NEB)	25 μl
i5-primer (25 μM)	2 μl
i7-primer (25 μM)	2 μl
Molecular biology water (46000-CM, Corning)	11 μl

Temperature profile

Step 1: 72°C for 5 min
Step 2: 98°C for 30 s
[Step 3: 98°C for 10 s
Step 4: 63°C for 30s
Step 5: 72°C for 1 min]
(Repeat steps 3-5 for a total number of 8 cycles)
Step 6: keep at 12°C.

9. Add 250 μl Buffer PB (Qiagen) and 10 μl Na-acetate (3 M, pH = 5.2). Purify using MinElute PCR Purification Kit (28004, Qiagen) and elute in 20 μl EB.
10. Size selection with Ampure XP beads (A63880, Beckman Coulter). Add 180 μl EB and 110 μl beads (0.55 x sample volume) and mix 10 times by pipetting.
11. Incubate 5 min at room temperature. Separate on magnetic stand.
12. Transfer supernatant (300 μl) to a new tube.
13. Add 190 μl beads (1.5 x sample volume).
14. Incubate 5 min at room temperature. Separate on magnetic stand.
15. Wash beads 2 times with 200 μl EtOH (70 %).
16. After the second wash resuspend beads in 20 μl EB by pipetting 10 times.
17. Separate on magnetic stand and transfer 17 μl of final library to a new tube.

Quantification of library and quality check

18. Quantify final libraries using Qubit (1 μ l/sample, Qubit dsDNA HS Assay Kit (Q32851), Thermo Fisher Scientific) and check for library size distribution using 4200 TapeStation (High Sensitivity D1000 ScreenTape (5067-5584) and Reagents (5067-5583), Agilent Technologies).
19. Optional: Quantify final library using the KAPA DNA Quantification Kit (KK4953, KAPA Biosystems)
20. Samples are now ready for sequencing on Illumina platform e.g. HiSeq4000 (50 bp PE, Nextera v2 libraries).

Appendix

A. PCR Primer

PCR primers were ordered from IDT. Primer sequences were the same as in the Nextera Index XT Kit v2 (FC-131-2001, Illumina).

i5-Primer

Name	Index sequence	Full sequence
S502	CTCTCTAT	AATGATACGGCGACCACCGAGATCTACACCTCTCTATTCGTCTCGGCAGCGTC
S503	TATCCTCT	AATGATACGGCGACCACCGAGATCTACACTATCCTCTTCGTCTCGGCAGCGTC
S505	GTAAGGAG	AATGATACGGCGACCACCGAGATCTACACGTAAGGAGTCGTCTCGGCAGCGTC
S506	ACTGCATA	AATGATACGGCGACCACCGAGATCTACACACTGCATATCGTCTCGGCAGCGTC
S507	AAGGAGTA	AATGATACGGCGACCACCGAGATCTACACAAGGAGTATCGTCTCGGCAGCGTC
S508	CTAAGCCT	AATGATACGGCGACCACCGAGATCTACACCTAAGCCTTCGTCTCGGCAGCGTC
S510	CGTCTAAT	AATGATACGGCGACCACCGAGATCTACACCGTCTAATTCGTCTCGGCAGCGTC
S511	TCTCTCCG	AATGATACGGCGACCACCGAGATCTACACTCTCTCCGTCGTCTCGGCAGCGTC

i7-Primer

Name	Index sequence	Full sequence
N701	TCGCCTTA	CAAGCAGAAGACGGCATAACGAGATTCGCCTTAGTCTCGTGGGCTCGG
N702	CTAGTACG	CAAGCAGAAGACGGCATAACGAGATCTAGTACGGTCTCGTGGGCTCGG
N703	TTCTGCCT	CAAGCAGAAGACGGCATAACGAGATTTCTGCCTGTCTCGTGGGCTCGG
N704	GCTCAGGA	CAAGCAGAAGACGGCATAACGAGATGCTCAGGAGTCTCGTGGGCTCGG
N705	AGGAGTCC	CAAGCAGAAGACGGCATAACGAGATAGGAGTCCGTCTCGTGGGCTCGG
N706	CATGCCTA	CAAGCAGAAGACGGCATAACGAGATCATGCCTAGTCTCGTGGGCTCGG
N707	GTAGAGAG	CAAGCAGAAGACGGCATAACGAGATGTAGAGAGGTCTCGTGGGCTCGG
N710	CAGCCTCG	CAAGCAGAAGACGGCATAACGAGATCAGCCTCGGTCTCGTGGGCTCGG
N711	TGCCTCTT	CAAGCAGAAGACGGCATAACGAGATTGCCTCTTGTCTCGTGGGCTCGG
N712	TCCTCTAC	CAAGCAGAAGACGGCATAACGAGATTCCTCTACGTCTCGTGGGCTCGG