

# High-Pressure Liquid Chromatography of Sialic Acids on a Pellicular Resin Anion-Exchange Column with Pulsed Amperometric Detection: A Comparison with Six Other Systems<sup>1</sup>

Adriana E. Manzi,<sup>2</sup> Sandra Diaz, and Ajit Varki<sup>3</sup>

*Division of Hematology-Oncology, Department of Medicine, The Cancer Biology Program, UCSD Cancer Center, and the San Diego Veterans Medical Center, University of California at San Diego, La Jolla, California 92093*

Received December 28, 1989

**A wide variety of different sialic acids have been reported in nature. Following their release and purification, detection and quantitation of these molecules is now possible by a number of techniques. We and others have previously reported high-pressure liquid chromatography separation of sialic acids with several different columns, elution methods, and detection techniques. We report here a new method for the separation of sialic acids at neutral pH on a Carbowac PA-1 anion-exchange column of pellicular resin, with pulsed amperometric detection following postcolumn addition of alkali. The major advantages of this system are the separation of a variety of sialic acids, sensitive detection (into the picomole range), and the relative ease of use for preparative purposes. Using a set of defined sialic acid standards, this method is compared and contrasted with six other HPLC methods previously described by us and by others. The advantages and disadvantages of each system are also addressed. In the final analysis, no single method is adequate to completely separate and quantitate all of the known sialic acids. However, used in appropriate combinations, these methods allow exploration of the biology of sialic acids in a manner heretofore not possible.** © 1990 Academic Press, Inc.

Press, Inc.

<sup>1</sup> This research was supported by RO1-GM32373 from the USPHS (to A.V.). A.M. was partially supported by a postdoctoral fellowship from the Consejo Nacional de Investigaciones Científicas y Técnicas, Argentina. A.V. is the recipient of a faculty research award from the American Cancer Society (FRA-295).

<sup>2</sup> Current address: Departamento de Química Orgánica, Facultad de Ciencias Exactas y Naturales, Universidad de Buenos Aires, Ciudad Universitaria, Pabellón 2, 1428 Buenos Aires, Argentina.

<sup>3</sup> To whom correspondence should be addressed at UCSD, Q-063, La Jolla, CA 92093.

The sialic acids are a family of 9-carbon carboxylated 2-keto sugars usually found as the terminal residues of eucaryotic oligosaccharides. More than 25 different kinds of modified sialic acids have now been reported in nature. The most commonly occurring sialic acid is *N*-acetylneuraminic acid (5-acetamido-3,5-dideoxy-D-glycero-D-galactononulosonic acid, or Neu5Ac,<sup>4</sup> which is also presumed to be the biosynthetic precursor of all the other molecules (1,2). One modification is the substitution of one of the hydrogens on the *N*-acetyl group by a hydroxyl group, giving rise to *N*-glycolylneuraminic acid (Neu5Gc). Another is the replacement of the amino group at the 5 position with a hydroxyl group, giving rise to 2-keto-3-deoxynonulosonic acid (KDN) (3). Most of the other sialic acids arise from O-substitution of one or more of the hydroxyl groups of Neu5Ac, Neu5Gc, or KDN with acetyl, methyl, lactyl, or sulfate groups. Unsaturated and dehydro forms of sialic acids have also been reported in nature (see Fig. 1, reviewed in Refs. (1,2,4). These modifications show tissue-specific and developmentally regulated distribution (5-9) and are known to affect a wide spectrum of biological phenomena (10-17). In order to further explore the biology of sialic acid modifications it is necessary to release and purify these compounds from biological sources in an intact state (18), and to identify and quantitate the various sialic acids present in small quantities in complex biological mixtures.

Characterization and quantitation of purified sialic acids can be achieved by thin-layer chromatography

<sup>4</sup> Abbreviations used: The various sialic acids are designated according to Schauer and others (1) using combinations of Neu (neuraminic), Ac (acetyl), Me (methyl), and Gc (glycolyl), e.g., Neu5,7Ac<sub>2</sub> is 7-*O*-acetyl-*N*-acetylneuraminic acid (see Table 1 for a complete list). Other abbreviations used: BSM, bovine submaxillary mucin; DMB, 1,2-di-

(19), gas-liquid chromatography (19), colorimetry (20), or colorimetry combined with HPLC (21). Structural analysis has been performed by gas-liquid chromatography-mass spectrometry (GLC-MS) (22,23) as well as by NMR spectroscopy (19,24), provided the necessary amount of material can be obtained. GLC and GLC-MS require derivatization of the purified sample prior to analysis (19,22) and NMR spectroscopy requires highly purified samples and is not suitable for the analysis of complex mixtures (24). We have recently also demonstrated quantitation and identification of sialic acids in partially purified, complex biological mixtures by fast atom bombardment-mass spectrometry (FAB-MS).<sup>5</sup>

In general, HPLC has the advantage of allowing rapid and direct quantitation of underivatized samples, characterization of a given sample (number of components), and tentative identification (by comparison with the elution times of known standards and/or coinjection with them). In some cases it is also possible to perform preparative chromatography under similar conditions. Different HPLC systems for the characterization of sialic acids have been reported (25-32). Most previously described approaches suffer from lack of specificity in the detection of sialic acids (absorbance at 200-210 nm), and the overlapping elution positions of several of the compounds. More specific detection requires prior derivatization (28,33).

Recently, high-performance anion-exchange chromatography with pellicular resins has been shown to separate a variety of positional isomers of neutral monosaccharides, oligosaccharides, and glycopeptides (34,35) as well as sialylated and phosphorylated oligosaccharides (36,37). Elution at high pH allows separation of the oligosaccharides as their oxyanions. This is coupled to pulsed amperometric detection (PAD), an electrochemical method that makes use of the ionization of hydroxyl groups at high pH (34,38), and is characterized by high sensitivity and relative specificity for compounds with hydroxyl groups. PAD enables detection of 10-100 pmol of each carbohydrate in a mixture, although relative detector responses could be different for each component. We have explored the use of high-performance anion-exchange chromatography on pellicular resins with pulsed amperometric detection for the analysis of a variety of sialic acids. To protect labile *O*-acetyl esters, separation was achieved with a pH gradient of acetate ion below neutral pH, and postcolumn addition of alkali was used to prepare the eluate for the PAD detector. In this study, we directly compare this new method with six other previously described HPLC systems and point out advantages and disadvantages of each.




FIG. 1. The sialic acids. The nine-carbon backbone common to all known sialic acids is shown in a stylized chair conformation. Substituents, linkages, and anhydro forms reported to date are indicated.

## MATERIALS AND METHODS

### *Sialic Acid Standards (see also Table 1)*

[4-<sup>14</sup>C]Neu5Ac and [9-<sup>3</sup>H]Neu5Ac were prepared from labeled [1-<sup>14</sup>C] and [6-<sup>3</sup>H]ManNAc (American Radiolabelled Chemicals), as previously described (39); [*O*-acetyl-<sup>3</sup>H]Neu5,7Ac<sub>2</sub> and Neu5,9Ac<sub>2</sub>, were prepared from rat liver Golgi vesicles labeled with [acetyl-<sup>3</sup>H]acetyl-coenzyme A (ICN) as described elsewhere (32); [acetyl-<sup>3</sup>H]- and [acetyl-<sup>14</sup>C]Neu5Ac were synthesized by treatment of neuraminic acid β-methyl glycoside (Neu2βMe) with a 1:1 (v/v) mixture of [2-<sup>3</sup>H]- or [1-<sup>14</sup>C]acetic anhydride (in toluene) and pyridine, for 24 h at RT. The resulting <sup>3</sup>H- or <sup>14</sup>C-labeled *N*-acetylneuraminic acid β-methyl glycosides (carrying some *O*-acetyl esters) were converted into *N*-acetylneuraminic acid by acid hydrolysis with 0.1 N HCl for 1 h at 80°C (which also removed *O*-acetyl esters), and purified by HPLC using System I. The neuraminic acid methyl ester, methyl-β-glycoside (Neu1,2Me<sub>2</sub>) was prepared by treating Neu2βMe in methanol with diazomethane in ether until a persistent faint yellow color was obtained, followed by evaporation with nitrogen. *N*-Acetylneuraminic acid (>99% purity) was from Kantoishi Pharmaceutical (MECT) Co., Tokyo, Japan, and *N*-glycolylneuraminic acid (Neu5Gc) from Sigma. 9-Acetamido-Neu5Ac was kindly provided by Professor Rienhard Brossmer, University of Heidelberg. 2,3-Dehydro, 2,6-anhydro-Neu5Ac (Neu2en5Ac) was from Boehringer-Mannheim. Sialic acids from bovine submaxillary mucin, equine serum, or the total lipid extract of the starfish *Pisaster brevispinus* were released, purified, and fractionated as previously described (32). Precautions to prevent loss and migration of *O*-acetyl esters were as previously described (18,32).

### Chemicals

1,2-Diamino-4,5-methylenedioxybenzene dihydrochloride (DMB) and 3-aminopropyl-triethoxysilane

amino-4,5-methylenedioxybenzene; ES, equine serum; KDN, 2-keto-3-deoxynonulosonic acid; PAD, pulsed amperometric detection; TBA, 2-thiobarbituric acid; RT, room temperature.

<sup>5</sup> Manzi, A., et al. (1990) *J. Biol. Chem.*, in press.

TABLE 1  
Sialic Acids Standards Used in This Study

Sialic acid	Abbreviation	Source	Structural confirmation
<b>Radioactive</b>			
[9- <sup>3</sup> H] <i>N</i> -Acetylneuraminic acid	[9- <sup>3</sup> H]Neu5Ac	Biosynthetic	Comigrates with authentic standard in several systems
[4- <sup>14</sup> C] <i>N</i> -Acetylneuraminic acid	[4- <sup>14</sup> C]Neu5Ac	Biosynthetic	As above
[acetyl- <sup>3</sup> H] <i>N</i> -Acetylneuraminic acid	[Ac- <sup>3</sup> H]Neu5Ac	Synthetic	As above
[acetyl-1- <sup>14</sup> C] <i>N</i> -Acetylneuraminic acid	[Ac-1- <sup>14</sup> C]Neu5Ac	Synthetic	As above
[9- <i>O</i> -acetyl- <sup>3</sup> H] <i>N</i> -Acetyl-9- <i>O</i> -acetylneuraminic acid	[9- <i>O</i> -Ac- <sup>3</sup> H]Neu5,9Ac <sub>2</sub>	Biosynthetic	As above
[7- <i>O</i> -acetyl- <sup>3</sup> H] <i>N</i> -Acetyl-7- <i>O</i> -acetylneuraminic acid	[7- <i>O</i> -Ac- <sup>3</sup> H]Neu5,7Ac <sub>2</sub>	Biosynthetic	As above
<b>Nonradioactive</b>			
<i>N</i> -Acetylneuraminic acid	Neu5Ac	Commercial	FAB-MS
<i>N</i> -Glycolylneuraminic acid	Neu5Gc	Commercial	FAB-MS
<i>N</i> -Acetyl-7- <i>O</i> -acetylneuraminic acid	Neu5,7Ac <sub>2</sub>	BSM	FAB-MS
<i>N</i> -Acetyl-9- <i>O</i> -acetylneuraminic acid	Neu5,9Ac <sub>2</sub>	BSM	FAB-MS
<i>N</i> -Acetyl-4- <i>O</i> -acetylneuraminic acid	Neu4,5Ac <sub>2</sub>	ES	FAB-MS
<i>N</i> -Acetyl-7(8),9-di- <i>O</i> -acetylneuraminic acid	Neu5,7(8),9Ac <sub>3</sub>	BSM	FAB-MS
<i>N</i> -Acetyl-7,8,9-tri- <i>O</i> -acetylneuraminic acid	Neu5,7,8,9Ac <sub>4</sub>	BSM	FAB-MS
<i>N</i> -Glycolyl-9- <i>O</i> -acetylneuraminic acid	Neu5Gc9Ac	BSM	FAB-MS
<i>N</i> -Glycolyl-4- <i>O</i> -acetylneuraminic acid	Neu4Ac5Gc	ESM	FAB-MS
<i>N</i> -Acetyl-2-deoxy-2,3-didehydroneuraminic acid	Neu2en5Ac	Commercial	FAB-MS
<i>N</i> -Glycolyl-8- <i>O</i> -methylneuraminic acid	Neu5Gc8Me	Starfish	FAB-MS
<i>N</i> -Glycolyl-8- <i>O</i> -methyl-9- <i>O</i> -acetylneuraminic acid	Neu5Gc8Me9Ac	Starfish	FAB-MS
<i>N</i> -Glycolyl-8- <i>O</i> -methyl-7,9-di- <i>O</i> -acetylneuraminic acid	Neu5Gc8Me7,9Ac <sub>2</sub>	Starfish	FAB-MS
Neuraminic acid $\beta$ -methylglycoside	Neu2Me	Synthetic	
Neuraminic acid methyl ester $\beta$ -methylglycoside	Neu1,2Me <sub>2</sub>	Synthetic	
<i>N</i> -Acetylneuraminic acid methyl ester $\beta$ -methylglycoside	Neu1,2Me <sub>2</sub> 5Ac	Synthetic	
<i>N</i> -Acetylneuraminic acid $\beta$ -methylglycoside	Neu2Me5Ac	Synthetic	
<i>N</i> -Acetyl-9- <i>O</i> -acetamidoneuraminic acid	Neu5Ac9NHAc	Synthetic	
<i>N</i> -Acetyl-4,8-anhydroneuraminic acid	Neu5Ac4,8anhydro	ESM	FAB-MS

(aminopropylsilane) were from Sigma. Diazomethane in ether was kindly provided by Dr. Alan Hoffman, Department of Medicine, University of California, San Diego. All other chemicals were of reagent grade and were purchased from commercial sources.

#### Analysis and De-*O*-acetylation of Free Sialic Acids

The 2-thiobarbituric acid (TBA) assay was used to measure free sialic acids (20). Analytical de-*O*-acetylation (18) to quantitate interference with the assay by *O*-acetyl substituents at the 7-, 8-, or 9-positions was carried out with 0.1 M NaOH for 30 min at 37°C, and quenched by addition of 0.2 M HCl (5  $\mu$ l for 10  $\mu$ l of reaction mixture) (32).

#### Induced Migration of *O*-Acetyl Esters

To induce migration of *O*-acetyl groups from the 7- to the 9-position with minimal loss of the esters, samples were treated with 1 ml of freshly prepared 10 mM NH<sub>4</sub>OH for 30 min at room temperature and dried down with a shaker evaporator (32).

#### Thin-Layer Chromatography

Sialic acids were analyzed by thin-layer chromatography on precoated cellulose plates using *n*-butanol:*n*-pro-

panol:0.1 M HCl, 1:2:1 (v/v/v) as developing solvent (40). The periodate/resorcinol reagent, as described by Jourdan *et al.* (41) was used for staining.

#### Derivatization with DMB

As described by Hara *et al.* (28), samples were heated at 50°C for 2.5 h in 7 mM DMB, 0.75 M  $\beta$ -mercaptoethanol, and 18 mM sodium hydrosulfite in 1.4 M acetic acid, in foil-covered tubes. The reagent is stable for up to 1 week at 4°C in the dark. Derivatized samples could be stored in the dark at 4°C for up to 72 h without excessive deterioration.

#### HPLC Columns

The columns used were from the following sources: CarboPac PA-1 (250  $\times$  4 mm i.d., and 250  $\times$  9 mm i.d.) anion-exchange columns and CarboPac PA-1 guard column, Dionex Corp.; Aminex HPX-72S (300  $\times$  7.8 mm i.d.; particle size, 11  $\mu$ m), Aminex HPX-87H (300  $\times$  7.8 mm; particle size, 9  $\mu$ m), and Aminex A-29 (40  $\times$  4.6 mm i.d.), Bio-Rad; Micropak AX-5 (300  $\times$  4 mm i.d.; particle size, 5  $\mu$ m), Varian; Versapak C-18 (250  $\times$  4.1 mm i.d.; particle size, 10  $\mu$ m), Alltech Assoc.; and TSK gel ODS-120-T (250  $\times$  4.6 mm i.d.; particle size, 5  $\mu$ m), TosoHaas.

### HPLC Pumps

A QIC reciprocating single-piston pump (Dionex Corp.) was used for the Carbowac column. All other analyses used Spectra-Physics 8700 and 8700XR pumps.

### Detectors

The detectors used were: SP-8440 uv-visible detector, Spectra Physics; FD-300 dual monochromator fluorescence detector (fitted with a 24  $\mu$ l flow-cell), Spectro Vision, Inc.; Ionchrom pulsed amperometric detector, Dionex Corp.; and Flow One Beta radioactive flow detector, Radiomatic Instruments.

### HPLC Solvents, Buffers, and Elution Procedures

Deionized water was passed through a five-stage Milli-Q Plus system, Millipore Corp. All solvents were HPLC grade, and all the buffers used for HPLC were prepared with Milli-Q purified water and filtered using a 0.45- $\mu$ m Nylon 66 membrane, Alltech Associates, Inc. Certified 50% (19.6 M) sodium hydroxide solution, with less than 0.05% sodium carbonate (Fisher Scientific) was used to prepare the 300 mM solution for postcolumn addition of base, prior to detection by the PAD.

#### System I—Ion-Exchange Chromatography on Pellicular Resin with PAD Detection

Carbowac PA-1 columns were eluted at RT in the isocratic mode for 5 min with 5 mM sodium acetate solution, followed by a 30-min linear gradient to 50% 5 mM sodium acetate, 50% 5 mM acetic acid. The column was then washed with 100% 5 mM acetic acid for 10 min, and with 100% 5 mM sodium acetate for a further 15 min. The flow was 1 ml/min for the analytical column (4  $\times$  250 mm) and 5 ml/min for the semipreparative column (9  $\times$  250 mm). Typical back pressure was 950–1400 psi. For analytical runs, the column effluent was mixed with 0.3–0.4 ml/min of 300 mM NaOH pumped under helium pressure, and monitored with the PAD using the following settings:  $E_1 = 0.05$ ;  $E_2 = 0.6$ ;  $E_3 = -0.8$ ;  $T_1 = 2$ ;  $T_2 = 2$ ;  $T_3 = 5$ ; response time = 1; range = 100 or 300 nA. In preparative runs, 0.25-min fractions were collected and checked by the TBA assay. Fractions were pooled accordingly and desalted with a Dowex 50 ( $H^+$ ) column, washed with water. Pooled fractions were recovered by lyophilization, and analyzed by TLC on cellulose plates.

#### System II—Amine Adsorption/Ion Suppression Chromatography

The Micropak AX-5 column was eluted in the isocratic mode with acetonitrile:water:0.25 M  $NaH_2PO_4$  (72:18:10), at 1 ml/min and RT (25). The proportions of the solvents depended upon the separation desired (see

Results). The effluent was monitored for uv absorption at 200–205 nm, or for radioactivity according to the sample. To detect radioactivity with maximum efficiency, it was necessary to use a high ratio of cocktail to effluent for the flow detector, or to evaporate the solvent from each collected fraction to be counted, and then add water and scintillation cocktail (1:10).

#### System III—Reverse-Phase Chromatography of Derivatized Sialic Acids

Sialic acids derivatized with DMB were applied to a TSK-ODS 120T column eluted isocratically with 0.9 ml/min of acetonitrile:methanol:water (9:7:84, v/v/v) at RT (28). Fluorescence of the eluate was monitored at an excitation wavelength of 373 nm and an emission wavelength of 448 nm; the high voltage was set at 900, the range at 500, and the response time at 0.5 s.

#### System IV—Anion-Exchange Chromatography on Aminex HPX-72S

The column was eluted isocratically with 0.1 M sodium sulfate solution (31,32). Depending on the particular column, the flow rate and working temperature were 1.0 ml/min at RT or 0.8 ml/min at 40°C (due to differences in the physical properties of the resin from batch to batch). The effluent was monitored for uv absorbance (200 nm) or for radioactivity. For preparative chromatography, fractions were collected beyond the uv detector directly into tubes on ice. Pooled fractions were diluted 20-fold in ice-cold distilled water and passed directly over a Dowex 3-X4A column (formate form) equilibrated in 10 mM sodium formate, pH 5.5, maintained at 4°C. Sialic acids were then eluted with 10 column vol of 1 M formic acid, and recovered by lyophilization.

#### System V—Anion-Exchange Chromatography on Aminex A-29

Elution was in the isocratic mode at RT at 0.5 ml/min flow, using 0.5–0.75 mM sodium sulfate (27,30). The exact concentration of sulfate depended upon the condition of the individual column. The effluent was monitored for uv absorbance (200 nm) or for radioactivity.

#### System VI—Strong Cation-Exchange Chromatography

The Aminex HPX-87H column was eluted isocratically with 0.006 M sulfuric acid at RT and 0.5 ml/min flow (29). The column effluent was monitored by means of a flow detector or by uv at 205 nm.

#### System VII—Reverse-Phase/Ion-Pair Chromatography

The eluting buffer was 0.4% tetrabutylammonium formate, pH 4.5, at 1.0 ml/min in the isocratic mode, at RT (25). The system was used only for radioactive samples,

and the radioactivity was continuously monitored by a flow detector, or by collecting fractions. For preparative runs, pooled fractions were passed directly over a 1 ml Dowex 50 ( $H^+$  form) column in water, washed through with 4 ml of water, and taken to dryness to remove the formic acid generated.


## RESULTS AND DISCUSSION

### *System I—CarboPac PA-1 Strong Anion-Exchange Pellicular Resin Column with Pulsed Amperometric Detection*

**Principles.** This type of anion-exchange chromatography on a polymeric, nonporous, pellicular resin was recently introduced for the analysis of neutral saccharides, which are either partially or completely ionized at high pH (12–14) (Ref. (38)). Separation of neutral carbohydrates is achieved by adjusting the eluant pH (Refs. (34,35)), and anionic oligosaccharides by sodium acetate/sodium hydroxide gradients (36,37). Such systems have separated polysaccharides and monosaccharides and have resolved structural isomers of oligosaccharides from biological samples. The only previous report of sialic acid fractionation is the separation of *N*-acetylneuraminic acid (retention time 4.0 min) and *N*-glycolylneuraminic acid (retention time 10 min) in the isocratic mode using 100 mM sodium hydroxide/150 mM sodium acetate as the eluent (Dionex handbook).

Pulsed amperometric detection utilizes a repeating sequence of three applied potentials, applied for specific durations (see Dionex Handbook for detailed principles). The total current which results is the sum of (a) the carbohydrate oxidation current, (b) the current due to the charging of the electrode surface, and (c) the current caused by the oxidation of the gold electrode. Therefore, the current is measured only during the application of  $E_1$  (after a delay of several tens to hundreds of milliseconds), thus improving the signal-to-noise ratio. PAD has been shown to be linear over more than four- and one-half orders of magnitude (38). Approximately 1% of the sample is oxidized, and molecules with similar size and structure are reported to have similar response factors (38). On the other hand, substituents on hydroxyl groups are expected to affect the sensitivity of detection. The detector works only at very high pH. Thus, when sialic acids are separated at neutral pH, it is necessary to add alkali to the effluent (300 mM sodium hydroxide at 0.3–0.4 ml/min) before it enters the PAD cell.

**Elution characteristics.** Since sialic acids have a  $pK_a$  of approximately 2.0 they are bound to the resin at neutral pH, and it is necessary to use an anion to elute them. Sodium hydroxide was avoided because *O*-acetyl esters could be saponified during chromatography. Samples were loaded in water to ensure binding and submitted



**FIG. 2.** Elution profile of standard sialic acids on a CarboPac PA-1 column with pulsed amperometric detection. A mixture of Neu5Ac, Neu5,9Ac<sub>2</sub>, Neu5,7(8),9Ac<sub>3</sub>, Neu4,5Ac<sub>2</sub>, and Neu5Gc (5 nmol of each standard) was studied. Neu5Ac and Neu5,9Ac<sub>2</sub> overlap closely, giving one broad peak.

to a pH gradient elution with dilute sodium acetate and acetic acid (see Materials and Methods). This permitted separation of different *O*-acetylated sialic acid species, whose binding to the resin is relatively weaker with higher numbers of substituents. Pure sialic standards and a variety of known mixtures from biological sources were analyzed by this method. Figure 2 shows an example with a mixture of known sialic acid standards. While some of the mono- and di-*O*-acetylated species can be separated from their parent molecule Neu5Ac, the relative detector responses vary significantly with the presence of different numbers of free hydroxyl groups (see below). Figure 3 shows the analysis of a mixture of sialic acids from equine serum (ES). After base hydrolysis, the two *O*-acetylated sialic acids disappeared as expected (Fig. 3B). In addition to the increase in the parent compounds, a broad peak appears corresponding to 4,8-anhydro compounds, a phenomenon peculiar to base hydrolysis of 4-*O*-acetylated sialic acids (2). With loss of two hydroxyl groups, the anhydro sialic acids also give a poorer detector response.

**Detection limits.** The reported limit of sensitivity of the PAD detector for carbohydrates is about 50 pmol (34). We found that accurate detection of *N*-acetylneuraminic acid can be done with as little as 200 pmol, using the above elution conditions (with the range setting of 100 nA, attenuation of the integrator at 1024, and a signal-to-noise ratio of 2).

**Relative detector responses.** As expected, the detector response of *N*-glycolylneuraminic acid, bearing one more hydroxyl group than *N*-acetylneuraminic acid, is 15% higher. In contrast, the response of the detector for




FIG. 3. Elution profile of a mixture of sialic acids from equine serum on a CarboPac PA-1 column with pulsed amperometric detection. Samples were injected onto the column untreated (A) or after base hydrolysis (B). The elution positions are as indicated.

Neu5,9Ac<sub>2</sub> is only 46%, and for Neu5,7(8),9Ac<sub>3</sub> 45% that of the response for Neu5Ac. Some of the relative detector response factors ( $\mu\text{V}\cdot\text{s}/\text{pmol}$ ) calculated for sialic acid standards from at least five different runs in different range and attenuation settings are: Neu5Ac, 30,500; Neu5,9Ac<sub>2</sub>, 14,500; Neu4,5Ac<sub>2</sub>, 14,100; Neu5Gc, 35,400; Neu5,7(8),9Ac<sub>3</sub>, 13,700; Neu5,7,8,9Ac<sub>4</sub>, 19,600; and Neu2en5Ac, 13,500. The detection limits are higher for substituted and anhydro sialic acids, because there are a lower number of ionizable hydroxyl groups. However, because of the possibility of partial de-O-acetylation between the point of addition of postcolumn alkali and the PAD detector, the detector response factors may vary somewhat in different laboratories.

**Relative elution times.** Under these conditions we were able to separate standard mixtures of Neu5Ac and Neu5Gc, their mono-, di-, and tri-O-acetylated derivatives, and Neu2en5Ac, and its mono-O-acetyl ester. However, some overlap occurred between various compounds (e.g., Neu5Ac and Neu5,9Ac<sub>2</sub>). Mixtures of purified sialic acids obtained from BSM and from the total lipid extract of the starfish *Pisaster brevispinus* whose compositions were previously determined by FAB-MS<sup>5</sup> were also tested. In the latter case, the Neu5Gc8Me was separated from its mono- and di-O-acetyl esters (not shown). Table 2 shows retention times of all the analyzed sialic acids (including some synthetic analogs) rel-

ative to Neu5Ac. This column can also be used for the analysis of labeled sialic acids, detected by means of an on-line flow detector, or by collecting and counting fractions (see Fig. 4 for an example of the latter). In this case, postcolumn addition of alkali is obviously not necessary.

**Use for preparative purposes.** For preparative runs the larger semipreparative column was used, and a maximum of 2  $\mu\text{mol}$  of total sialic acids could be injected without overloading the column. PAD could not be used since the postcolumn sodium hydroxide would hydrolyze the O-acetyl esters. We therefore first checked the profile of the mixture to be separated using the analytical column, and the PAD. The semipreparative column was then checked and calibrated with labeled standards (a mixture of [4-<sup>14</sup>C]Neu5Ac, [<sup>3</sup>H-O-acetyl]Neu5,7Ac<sub>2</sub>, and Neu5,9Ac<sub>2</sub>) using the flow detector to monitor the profile (as in Fig. 4). Finally, the preparative chromatography was done by collecting fractions, and monitoring the profile with the TBA assay for sialic acids (with analytical de-O-acetylation). Different areas were pooled,

TABLE 2  
Relative Elution Times of Sialic Acids  
in the Different Systems

Sialic acid	<i>R</i> <sub>Neu5Ac</sub> in system						
	I	II	III	IV	V	VI	VII
Neu5Ac	1.00	1.00	1.00	1.00	1.00	1.00	1.00
Neu5Gc	1.17	1.50	0.84	1.20	1.33	0.93	1.00
Neu5,7Ac <sub>2</sub>	0.74	0.36	1.06	0.95	0.87	0.97	1.71
Neu5,9Ac <sub>2</sub>	0.95	0.35	1.57	1.30	1.47	0.97	1.71
Neu4,5Ac <sub>2</sub>	0.76	0.39	1.68	1.22	1.41		
Neu5,7(8),9Ac <sub>3</sub>	0.74	0.23	1.90	1.17	1.76		
Neu5,7,8,9Ac <sub>4</sub>	0.62	<sup>a</sup>	1.98	1.62	2.05		
Neu5Gc9Ac	1.06	0.49	0.88	1.57	1.70		
Neu4Ac5Gc	0.86	0.54	1.53	1.59	1.69		
Neu2en5Ac	2.21	0.66	<sup>b</sup>	1.54	1.68		
Neu5Gc8Me	0.82	0.75	1.03	0.94	0.89		
Neu5Gc8Me9Ac	0.67	0.36	1.70	1.21			
Neu5Gc8Me7,9Ac <sub>2</sub>	0.59	0.29	2.20	1.47			
Neu2Me	0.17	0.13	<sup>b</sup>	0.23			
Neu1,2Me <sub>2</sub>	0.40	0.85	<sup>b</sup>	0.46			
Neu2Me5Ac	0.73	1.02	<sup>b</sup>	0.34	0.73		
Neu5Ac9NHAc	0.85	0.68		0.96			
Neu5Ac4,8anhydro	2.15	0.58	1.39	2.04			
Acetic acid	0.11	0.08			1.70		

*Note.* Values in roman type are derived from our studies. Values in italic are from others reported in the literature. Typical *R*<sub>Neu5Ac</sub> values for the various systems were as follows: System I, 18–24 min; System II(72:18:6), 21–24 min; System III, 12–13 min; System IV, 15–17 min; System V, 5–6 min; System VI, 11–12 min; System VII, 12–13 min. Note that relative elution times can vary somewhat with specific columns and temperatures, particularly with the isocratic systems. However, the order of elution of the various compounds relative to each other and to Neu5Ac is always maintained.

<sup>a</sup> Does not react with DMB.

<sup>b</sup> Elutes close to salt peak.




FIG. 4. Elution profile of labeled sialic acids on the CarboPac PA-1 column. A mixture of [*O*-acetyl- $^3\text{H}$ ]Neu5,7Ac $_2$  and [*O*-acetyl- $^{14}\text{C}$ ]Neu5,9Ac $_2$  with an internal standard of [ $^{14}\text{C}$ ]Neu5Ac was injected on the analytical CarboPac PA-1 column, and 0.5-min fractions were collected and counted.

passed over a Dowex 50 ( $\text{H}^+$ ) column (1 ml in water) to eliminate the acetate, lyophilized, and monitored by TLC on cellulose plates. One example of such preparative chromatography is shown in Fig. 5, for the mixture of sialic acids obtained from BSM. The recovery of total sialic acids was 87% in this case. As shown by the TLC (see Fig. 5B) there are four general areas in the profile: (i) di- (and probably tri-) *O*-acetylated sialic acids area, (ii) mono-*O*-acetylated Neu5Ac area, (iii) Neu5Ac plus mono-*O*-acetylated Neu5Gc area, and (iv) Neu5Gc area. Although several of the pools obtained contain more than one component, all showed substantial enrichment in one class of molecules. Thus, it is possible to do several preparative runs with 2  $\mu\text{mol}$  of the starting mixture, pool desired areas from each run, and rerun the combined pools for further purification.

**Advantages.** A major advantage of this system is the preparative separation of a variety of substituted sialic acids without exposure to conditions that cause loss or migration of *O*-acetyl esters. The acetate can be easily removed by Dowex 50 ( $\text{H}^+$  form) followed by lyophilization, again without effects upon the *O*-acetyl esters. Total recovery of sialic acids is also quite good (ca. 90%). Another advantage is analytical separation with detection at a high sensitivity (minimum of 200 pmol for Neu5Ac). Finally, this is one of the systems capable of separating the 7- and 9-mono-*O*-acetyl isomers of Neu5Ac and Neu5Gc.

**Limitations and disadvantages.** For preparative purposes, it is not possible to use the PAD to monitor the profile since the *O*-acetylated sialic acids are destroyed by the high concentration of sodium hydroxide. Thus, it is necessary to use a more laborious procedure (TBA reaction followed by TLC analysis) to monitor the preparation. The capacity of the preparative column currently marketed is somewhat limited (maximum 2  $\mu\text{mol}$

of sialic acids). Also, the relative detector responses of the sialic acids in the PAD depend upon the number of free hydroxyl groups. The possibility of variable de-*O*-acetylation at the point of addition of postcolumn alkali (prior to entry into the PAD detector) further complicates the determination of accurate detector responses.

**Column maintenance.** If the cleaning procedure with 100% of 5 mM acetic acid is not used in between runs the capacity of the column is significantly diminished after about 60 injections of 10 nmol each. When complete regeneration of the resin becomes necessary, the following procedure is useful: (i) run 3 column vols of 1 M HCl, (ii) wash well with water, (iii) run 3 column vols of 1 M NaOH, and (iv) wash well with water. After this procedure the retention times return to their normal values. However, with some columns the decreased resolution is not completely recovered.

#### Comparison of the Dionex System with other HPLC Systems

The set of sialic acid standards described above were also analyzed using six other previously described HPLC systems for sialic acids. In several cases modifications in described procedures were also devised to improve anal-




FIG. 5. Elution profile of sialic acids from bovine submaxillary mucin on a semipreparative CarboPac PA-1 column. Sialic acids from bovine submaxillary mucin (2  $\mu\text{mol}$ ) were chromatographed and fractions (0.25 min) collected. (A) The elution profile determined by the TBA reaction with and without prior de-*O*-acetylation. (B) The TLC profile of the pooled fractions on cellulose plates, using equal amounts (10 nmol) of sialic acids from each pool. Standards of Neu5Ac (extreme left lane) and Neu5Gc (extreme right lane) were also included. The plates were developed with resorcinol reagent. Because of poor photographic reproduction, the plates were photocopied and scanned using an Apple scanner. The location of the individual bands is thus indicated accurately.

ysis. The results for each system are presented below along with an outline of basic principles, advantages, and disadvantages in each case.

*System II—Amine Adsorption/Ion Suppression  
Chromatography: Micropak AX-5 Column  
with UV Detection*

This approach is based upon hydrogen bonding between the hydroxyl groups of saccharides and the amine functions of the stationary phase (42,43). In order to fractionate anionic molecules, an ion pair is added to the mobile phase to suppress the ionic effects, while retaining hydrogen bonding (44). We have previously reported use of this type of chromatography with a Micropak AX-5 column, for the analysis of some sialic acids (25). In this study, we report further experience with this system, and the retention times of several additional sialic acids. We obtained good separations and sharp peaks in the isocratic mode with a mixture of acetonitrile:water:0.25 M monobasic sodium phosphate at 1 ml/min. The relative ratio of the solvents can be varied depending on the kind of sialic acid mixture to be separated. It is necessary to maintain a minimum concentration of 10% phosphate and a maximum working percentage of acetonitrile (72%) above which the phosphate starts precipitating. A 64:26:10 (v/v/v) ratio of acetonitrile:water:0.25 M sodium phosphate is useful for the separation of Neu5Ac and Neu5Gc. However, for O-acetylated sialic acids it is necessary to increase the percentage of acetonitrile to increase their retention times (72:18:10, v/v/v). Chromatography is monitored by absorption at 200 nm, which requires extreme purity of the samples and reagents. Residual amounts of purification reagents such as acetic or formic acids produce an absorption peak close to the void volume. Typical profiles on this column for some radioactive sialic acids have previously been published (25). Table 2 shows the relative retention times for several sialic acids previously reported by us and those for some sialic acids and sialic acid analogs not previously reported. The lower limit of confident detection was 2 nmol, and the molar responses are similar for all the sialic acids.

The same analytical column can be used for preparative purposes (up to 400 nmol of sialic acids with a good resolution of the peaks). The eluant was monitored by uv and the peaks were collected directly on ice, diluted 10-fold with cold water, and purified by ion-exchange chromatography on Dowex 3-X4A (formate) at 4°C. The column was washed with 10 vol of 10 mM formic acid and eluted with 7 vol of 1 M formic acid, and sialic acids were recovered by lyophilization. We obtained ~65% recovery of sialic acids, but some loss and migration of O-acetyl groups occurred.

The various advantages and disadvantages of this system are summarized in Table 3. These columns also have

many other uses in the study of oligosaccharides (42,45,46), making them economically attractive. When a given analysis takes more than one day, it is recommended that a slow flow of the working buffer be maintained overnight. When the column loses resolution and retention is reduced, it can be cleaned by running 0.5 M phosphoric acid for at least 1 h and washing extensively with water before changing to the working buffer. More complete regeneration of the column can be achieved by injecting 3-aminopropyltriethoxysilane ( $3 \times 1$  ml) while running the column in 100% acetonitrile, washing extensively with the same solvent, and then with water and 0.5 M phosphate before use. This treatment restores the amino groups on the resin. The performance of the column is fully recovered after this treatment; however, the retention times may not be identical.

*Other similar systems.* Other amine-adsorption/ion-suppression based methods have been reported to be useful for the separation of sialic acids and sialic acid analogs, including Spherisorb-NH<sub>2</sub> columns (47,48). We have also explored use of a Varian Micropak-SP oligonucleotide column ( $300 \times 4$  mm), that can be run basically under the same conditions discussed for the Micropak AX-5, with similar results (not shown). However, this size of oligonucleotide column is no longer available.

*System III—Reverse-Phase Chromatography of  
Derivatized Sialic Acids:TSK-ODS 120T  
with Fluorometric Detection*

Hara *et al.* (33) have developed an HPLC separation of fluorescent derivatives of Neu5Ac and Neu5Gc, obtained by reaction with DMB in sulfuric acid and separated by reverse-phase chromatography on a Radial-Pak C18 cartridge with a mixture of methanol-acetonitrile-water. More recently, these authors formed the fluorescent derivatives in acetic acid at lower temperatures, completely avoiding migration and de-O-acetylation (28). An extensive list of other DMB-reacting substances that could be present in biological samples was also examined, all of which gave very weak responses and/or very long retention times under the recommended conditions (28). We have confirmed the value of this technique and extended it to the analysis of several additional sialic acid molecules (see Table 2). Figure 6 shows the elution profile of a complex mixture of sialic acids from BSM, with or without prior treatment with alkali. It can be seen that after induced migration with a very mild base, the peaks corresponding to the 7-O-acetylated molecules disappear because of migration of the esters to the 9-position (18,32,49). After more extensive base hydrolysis, there is almost complete elimination of the O-acetylated molecules, with an increase in the amount of the parent compounds, Neu5Ac and Neu5Gc.

The previously reported limits of detection were in the range of 100–200 fmol at a signal-to-noise ratio of 3 (28).



TABLE 3  
Comparison of Some Features of HPLC Systems for Sialic Acids

Feature	System						
	I	II	III <sup>a</sup>	IV	V	VI	VII
1. Sensitivity	Good	Poor	Best	Poor	Poor	Poor	NA <sup>b</sup>
2. Specificity of detection	Good	Poor	Best	Poor	Poor	Poor	NA <sup>b</sup>
3. Separation by number of hydroxyls or substituents	Good	Good	Good	Good	Good	Poor	Good
4. Separation of isomers	Good	Poor	Good	Good	Good	No	No
5. Preparative use	Best	Fair	NA	Good	Poor	NT	Good <sup>b</sup>
6. Avoids ester migration during purification	Good	Poor	Good	Poor	NT	NT	NT
7. Avoids ester loss during purification	Good	Fair	N/A	Fair	NT	NT	Good

Note. NA, not applicable, NT, not tested.

<sup>a</sup> Derivatized samples only.

<sup>b</sup> Radioactive samples only.

In our hands, confident detection and quantitation was possible with 2.5 pmol of any sialic acid. The original study reported no significant differences in detector responses between the various sialic acids (28). However,

we have some observed differences in the area of the peaks corresponding to the derivatives of different sialic acids, and found that it is probably due to the different degree of derivatization that each individual molecule reaches. In the prior study, the fluorescent peak profiles suggested that complete derivatization of the sialic acids might be occurring. However, analysis of labeled sialic acids by this method showed that the derivatization reaction does not go to completion. Thus, as shown in Fig. 7, different radioactive sialic acids gave different degrees of incorporation of radioactivity into the corresponding chromophore peak. A substantial portion of the radioactivity ran in other areas of the chromatogram, not corresponding to any fluorescence (note that the sample analyzed was identical to that shown in Fig. 4, and did not




FIG. 6. Elution profile of DMB-derivatized sialic acids on a TSK-ODS 120T column with fluorometric detection. A mixture of purified sialic acids from bovine submaxillary mucin (16.6 nmol) was derivatized with DMB after no treatment (A), induced migration of *O*-acetyl groups (B), or de-*O*-acetylation with base (C). The position of elution of known DMB-derivatives of sialic acids is indicated.




FIG. 7. Elution profile of DMB-derivatized radioactive sialic acids on a TSK-ODS 120T column. A mixture of [*O*-acetyl-<sup>3</sup>H]Neu5,7Ac<sub>2</sub> and Neu5,9Ac<sub>2</sub> with an internal standard of [<sup>14</sup>C]Neu5Ac (same sample as in Fig. 4) was derivatized with DMB and run on the TSK-ODS 120T column. Fractions (0.5 min) were collected and counted. The positions of elution of known fluorometric standards are indicated. The unknown radioactive peaks are indicated by a question mark.

contain any intrinsic labeled contaminants). Attempts to improve derivatization by increasing the reaction time to 4.5 h, or the concentration of reagent to 21 mM, did not result in an increase of the derivative peak, but of the third radioactive peak, indicating breakdown of the product increased (data not shown). However, as long as appropriate standards are run in parallel, the incomplete derivatization does not pose a practical problem in quantitation of unlabeled sialic acids.

The advantages and disadvantages of this system are summarized in Table 3. Preparative use of this column is not possible, because the samples are irreversibly derivatized. However, derivatized samples can be pooled, lyophilized, and rechromatographed after base treatment, to confirm O-acetylation (28). The original authors reported that the column could be used for more than 1000 injections (10  $\mu$ l each) if it was washed with methanol:water (1:1, v/v) at 1.0 ml/min for 30 min at the end of each day (28). However, in our hands, the pressure goes over the manufacturers recommended limit; we therefore wash at 0.5 ml/min. We also found it important to keep the amount of DMB reagent in each injection as low as possible to maintain column performance. Therefore, when very small amounts of material that need to be injected completely into the column are analyzed, we suggest a reduction in the total volume of the derivatization reaction.

#### *System IV—Anion-Exchange Chromatography on an Aminex HPX 72-S Column with uv Detection*

This column is normally used to separate weak organic bases or acids by a combination of ion-exclusion, ion-exchange, and partition chromatography on an 8% crosslinked quaternized polystyrene divinylbenzene copolymer in the sulfate form. When moderately concentrated sulfate solutions were used as the mobile phase, we found that this column separated sialic acids on the basis of the negative charge and the total number of hydroxyl groups (32). However, the presence of an extra hydroxyl group in the *N*-acetyl group (as in Neu5Gc) increased retention time relative to Neu5Ac (31,32). We have previously reported elution of sialic acids in the isocratic mode with 100 mM sodium sulfate, at room temperature and a flow rate of 1.0 ml/min (32). However, the manufacturer has recently reported a change in the physical characteristics of the resin, causing a need to change operating parameters. It is now necessary to increase the temperature (at least up to 40°C; the maximum temperature allowed by the resin is 70°C) and to reduce the flow rate (at least down to 0.8 ml/min) to achieve the same separation and to avoid compaction of the resin. About 2 nmol of each sialic acid can be confidently detected. As commented before for other HPLC methods, the relative detector responses can be assumed to be the same for the different sialic acid derivatives.

We previously reported the use of an Aminex HPX-72S column to separate Neu5,7Ac<sub>2</sub> and Neu5,9Ac<sub>2</sub> (32) and recently the utilization of the system for the analysis and preparative purification of complex mixtures of sialic acids from biological sources.<sup>5</sup> We report here also the retention times of a variety of other sialic acids and sialic acid analogs (see Table 2).

For preparative separations, the load limit is about 500 nmol of total sialic acids. Peaks are collected on ice directly from the effluent, immediately diluted with water (20-fold) and submitted to anion-exchange chromatography on Dowex 3-X4A (formate form) at 4°C. Sialic acids are eluted with 1 M formic acid and recovered by lyophilization. This procedure causes a low but variable loss of O-acetylated species because of migration and de-O-acetylation during post-HPLC purification. We have recently employed this method to completely separate pairs of sialic acids obtained by cellulose chromatography of total sialic acids from bovine submaxillary mucin.<sup>2</sup>

The major advantages and disadvantages of this method are summarized in Table 3. The retention time of acetate is quite long in this system (approx 30 min), eluting beyond all the sialic acids. This can allow the characterization of the sample and the quantitation of the percentage of O-acetyl esters present in the starting material, following base treatment. A major current limitation of this method is the variation in availability of a resin of the required quality to resist the high pressures involved in the analysis. While the previously available columns showed stable performance for over 250 injections, the physical properties of the packing material are no longer as consistent. It has become necessary to work at a lower flow rate and higher temperatures in order to avoid high back pressure. The higher pressure limit is set at 1000 psi because whenever the working back pressure exceeds this value, the subsequent increase in back pressure becomes exponential. When the limit is reached, it is possible to further increase the temperature (up to the limit of 70°C) or reduce the flow rate, and to recheck the retention times and resolution of standards in the new conditions. In order to extend the lifetime of the column it is also necessary to use a guard column. It is advisable to start heating the column with solvent flowing through it, but not to initially operate the column at a flow rate greater than 0.2 ml/min at ambient temperature. During the shutdown procedure the heater is first turned off and the flow rate is reduced to 0.2 ml/min and continued until the column reaches ambient temperature. The column can be cleaned with 5 or 20% acetonitrile in 0.1 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> and regenerated with 0.5 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>.

#### *System V—Anion-Exchange Chromatography on Aminex A-28 or A-29 with uv Detection*

Bio-Rad Aminex ion-exchange resins for HPLC are finely sized spherical beads of polystyrene divinylben-

zene copolymer with attached functional groups. Shukla and Schauer reported anion-exchange HPLC of sialic acids carried out on quaternary amine resin types, with 8% crosslinkage, and a particle size of  $11 \pm 2 \mu\text{m}$  (Aminex A-28) or  $5\text{--}8 \mu\text{m}$  (Aminex A-29) (27,30,50,51). These very short columns ( $40 \times 4.6 \text{ mm}$ ) are eluted in the isocratic mode with dilute sodium sulfate. Schauer and co-workers usually use 0.75 mM sodium sulfate at a flow rate of 0.5 ml/min and ambient temperature. In our hands, the resolution is improved when 0.50 mM sodium sulfate is used as the eluant at that flow rate. This difference might depend on the particular column that is used, because of the extremely small column volume. Each run takes approximately 10 min and can be monitored by on-line reading of the absorbance at 200 nm, or by monitoring of radioactivity. In our hands, the detection of 2 nmol of Neu5Ac is possible. However, Schauer's group reported lower limits of detection (27,50). The relative detector response can be considered the same for all compounds in the family that have an *N*-acyl group. Table 2 contains all the reported retention times for sialic acids and sialic acid analogs including our data and the data of Schauer.

The general advantages and disadvantages of this system are summarized in Table 3. The overlapping of peaks in complex mixtures is too high and the capacity of the column too low to permit preparative chromatography. A major advantage of this system is the short time required for the analysis, making it specially suitable for the study of enzyme reactions, where one aliquot of the enzyme assay mixture can be directly injected into the column. This usually gives information about both substrates and reaction products in one chromatographic run. Several such examples have been reported (27,30,52). The fast elution of substituted sialic acids also helps to avoid migration of acetyl groups and/or de-O-acetylation. However, because of the short time of the run, elution times are very close and it is easy to miss the presence of a minor component. Therefore, it is not recommended for the analysis of complex mixtures from biological sources.

The column can be used for the analysis of total enzyme reaction mixtures if small volumes are injected, but when crude enzymatic preparations are used it is advisable to previously purify the sialic acids. The column can be cleaned by running it backward for 30 min in 5% acetonitrile in the working buffer at 0.1 ml/min; the percentage of acetonitrile can be increased up to 30% if necessary. More complete regeneration is done with 0.2 M ammonium or sodium sulfate in water at 0.1 ml/min, at room temperature, overnight. However, in our hands the life of the column is relatively limited.

#### *System VI—Strong Cation-Exchange Chromatography: Bio-Rad HPX-87H Column with uv Detection*

The Aminex HPX-87H column is packed with cation-exchange resin in the hydrogen form. It was previously

used for the analysis and quantitation of Neu5Ac in serum of melanoma and breast carcinoma patients and was able to separate Neu5Ac and Neu5Gc (29). A mobile phase of 0.006 N sulfuric acid was previously used at a flow rate of 0.65 ml/min and  $42^\circ\text{C}$  for the separation of Neu5Ac, and Neu5Gc, with the uv detector set at 206 nm (29). In order to maintain the back pressure of the column within the permitted limits, we used a flow rate of 0.5 ml/min, at room temperature. The maximum reported sensitivity for Neu5Ac in this system is 0.8 nmol/ml (total injection volume =  $50 \mu\text{l}$ ) with a signal-to-noise ratio  $>3$  (29). Previously reported elution times are included in Table 2, together with some of our recent data in this system. Analysis of sialic acids on this column is limited only to Neu5Ac and Neu5Gc, since substituted sialic acids coelute with their parent molecules and may be destroyed by sulfuric acid. No major advantage of this system has been found (see Table 3). It provides a faster comparative analysis of sialic acids levels in serum hydrolysates than the TBA-HPLC assay (21) (because of the lack of the derivatization step), but is not better in this regard than the Aminex A-28 or A-29 system. The use for preparative purposes has many obvious limitations. Similar precautions to those for the HPX 72-S column are necessary to control the back pressure and to prevent compaction.

#### *System VII—Reverse-Phase/Ion-Pair Chromatography: C18 Column with Radioactive Detection*

Reverse-phase, ion-pair, column partition chromatography neutralizes the negatively charged carboxyl groups of the sialic acids with a positively charged lipophilic counterion (e.g., tetrabutylammonium) dissolved in the mobile phase, to form an ion pair that is able to partition into the lipophilic stationary phase of the column. The elution of the ion pair is produced by increasing the solvent strength (e.g., methanol concentration) of the mobile phase. We previously reported the use of a  $\mu\text{Bondapak C18}$  column (Waters) with TBA formate as a counterion for the analysis of free sialic acids (25). This system cannot be monitored by uv because of the high absorbance of the buffer, limiting analysis to radioactive samples. The relative elution times of the labeled sialic acids that have been analyzed by this method are included in Table 2. The system cannot separate Neu5Ac and Neu5Gc, but allows separation according to the number of *O*-acetyl groups, due to the increase in hydrophobicity.

Although the use of this system for preparative purposes has not been explored extensively, it is possible to easily eliminate the tetrabutylammonium formate by passing the eluant through a Dowex 50 ( $\text{H}^+$  form) column. At the present time, we see no obvious advantage to this method over other previously described methods. However, because it can only be used for radioactive si-

alic acids, we have not been able to examine a large number of standards. As with many other reverse-phase separations, it shows significant variability in elution time if the column is not thoroughly cleaned with organic solvents such as methanol and well-equilibrated with the ion-pair containing solution.

## CONCLUSIONS

In this study, we have reported a new HPLC method for the separation and analysis of sialic acids: anion-exchange chromatography on a pellicular resin column with pulsed amperometric detection. The high sensitivity of this type of detection gives the system the lowest detection limit available for the HPLC analysis of underivatized sialic acids (minimum of 200 pmol). The CarboPac PA-1 column was especially useful for preparative purposes, allowing the separation of a variety of substituted sialic acids without loss or migration of *O*-acetyl esters and with excellent recovery. Using a variety of natural and synthetic sialic acids, we have compared this method with all of the previously described HPLC methods for the analysis of these compounds. In four of these systems, we have made a direct comparison with the present method, using the same set of standards, and have reported new information to build upon the existing literature. Two other systems of lesser utility that we have experience with are mentioned primarily for completeness. We have mentioned the advantages and disadvantages of each system, and made general suggestions for their use (see Table 3 for a summary). We have also described the basic principles of each method and the precautions to be taken in their use. It should be noted that less than half of all the possible combinations of sialic acid modifications have been currently reported in the literature. Thus, it is possible that the use of these techniques will result in the discovery of new sialic acids, with different combinations of currently known substitutions, or new types of substitutions.

In general, the effect of substitutions on retention times is different for different systems. For example, *O*-acetyl esters retard elution in the HPX 72-S system, whereas free hydroxyl groups cause retardation in the AX-5 system. Thus, it is possible to roughly predict where a particular type of sialic acid might elute in the various systems. However, variations between isomers occur, making exact predictions impossible. It is also important to note that in all isocratic separations, differences in the exact retention times can be expected between different columns or different laboratory conditions. However, the order of elution of the different sialic acids is very consistent. In the final analysis, no single system is ideal for all aspects of the preparation and analysis of the various sialic acids. However, taken together, the various methods allow the detailed study of the wide variety of sialic acids that are found in nature.

## REFERENCES

1. Rosenberg, A., and Schengrund, C. (1976) in *Biological Roles of Sialic Acid* (Rosenberg, A., and Schengrund, C., Eds.), pp. 295-359, Plenum, New York/London.
2. Schauer, R. (1982) *Sialic Acids: Chemistry, Metabolism and Function*, Cell Biology Monographs, Vol. 10, Springer-Verlag, New York.
3. Nadano, D., Iwasaki, M., Endo, S., Kitajima, K., Inoue, S., and Inoue, Y. (1986) *J. Biol. Chem.* **261**, 11,550-11,557.
4. Schauer, R. (1982) *Adv. Carbohydr. Chem. Biochem.* **40**, 131-234.
5. Muchmore, E., Varki, N., Fukuda, M., and Varki, A. (1987) *FASEB J.* **1**, 229-235.
6. Constantine-Paton, M., Blum, A. S., Mendez-Otero, R., and Barnstable, C. J. (1986) *Nature (London)* **324**, 459-462.
7. Sparrow, J. R., and Barnstable, C. J. (1988) *J. Neurosci.* **21**, 398-409.
8. Levine, J. M., Beasley, L., and Stallcup, W. B. (1986) *Dev. Brain Res.* **27**, 211-222.
9. Schlosshauer, B., Blum, A. S., Mendez-Otero, R., Barnstable, C. J., and Constantine-Paton, M. (1988) *J. Neurosci.* **8**, 580-592.
10. Varki, A., and Kornfeld, S. (1980) *J. Exp. Med.* **152**, 532-544.
11. Higa, H. H., Rogers, G. N., and Paulson, J. C. (1985) *Virology* **144**, 279-282.
12. Rogers, G. N., Herrler, G., Paulson, J. C., and Klenk, H. D. (1986) *J. Biol. Chem.* **261**, 5947-5951.
13. Muchmore, E., and Varki, A. (1987) *Science* **236**, 1293-1295.
14. Cheresch, D. A., Reisfeld, R. A., and Varki, A. (1984) *Science* **225**, 844-846.
15. Higashi, H., Hirabayashi, Y., Fukui, Y., Naiki, M., Matsumoto, M., Ueda, S., and Kato, S. (1985) *Cancer Res.* **45**, 3796-3802.
16. Orskov, F., Orskov, I., Sutton, A., Schneerson, R., Lin, W., Egan, W., Hoff, G. E., and Robbins, J. B. (1979) *J. Exp. Med.* **149**, 669-685.
17. Corfield, A. P., Sander Wewer, M., Veh, R. W., Wember, M., and Schauer, R. (1986) *Biol. Chem. Hoppe Seyler* **367**, 433-439.
18. Varki, A., and Diaz, S. (1984) *Anal. Biochem.* **137**, 236-247.
19. Schauer, R. (1987) in *Methods in Enzymology* (Ginsburg, V., Ed.), Vol. 138, pp. 132-161, Academic Press, San Diego.
20. Warren, L. (1959) *J. Biol. Chem.* **234**, 1971-1975.
21. Powell, L. D., and Hart, G. W. (1986) *Anal. Biochem.* **157**, 179-185.
22. Casals Stenzel, J., Buscher, H. P., and Schauer, R. (1975) *Anal. Biochem.* **65**, 507-524.
23. Reuter, G., Pfeil, R., Stoll, S., Schauer, R., Kamerling, J. P., Ver-sluis, C., and Vliegthart, J. F. (1983) *Eur. J. Biochem.* **134**, 139-143.
24. Haverkamp, J., van Halbeek, H., Dorland, L., Vliegthart, J. F., Pfeil, R., and Schauer, R. (1982) *Eur. J. Biochem.* **122**, 305-311.
25. Diaz, S., and Varki, A. (1985) *Anal. Biochem.* **150**, 32-46.
26. Shukla, A. K., Scholz, N., Reimerdes, E. H., and Schauer, R. (1982) *Anal. Biochem.* **123**, 78-82.
27. Shukla, A. K., and Schauer, R. (1986) *Anal. Biochem.* **158**, 158-164.
28. Hara, S., Yamaguchi, M., Takemori, Y., Furuhashi, K., Ogura, H., and Nakamura, M. (1989) *Anal. Biochem.* **179**, 162-166.
29. Silver, H. K. B., Karim, K. A., Gray, M. J., and Salinas, F. A. (1981) *J. Chromatogr.* **224**, 381-388.
30. Schauer, R., Schroder, C., and Shukla, A. K. (1984) *Adv. Exp. Med. Biol.* **174**, 75-86.

31. Holzhauser, R., Faillard, H., Klose, W., Huber, W., Stickl, H., and Landthaler, M. (1988) *Klin. Wochenschr.* **66**, 540–544.
32. Diaz, S., Higa, H. H., Hayes, B. K., and Varki, A. (1989) *J. Biol. Chem.* **264**, 19,416–19,426.
33. Hara, S., Yamaguchi, M., Takemori, Y., Nakamura, M., and Ohkura, Y. (1986) *J. Chromatogr.* **377**, 111–119.
34. Hardy, M. R., Townsend, R. R., and Lee, Y. C. (1988) *Anal. Biochem.* **170**, 54–62.
35. Hardy, M. R., and Townsend, R. R. (1989) *Carbohydr. Res.* **188**, 1–7.
36. Townsend, R. R., Hardy, M. R., Hindsgaul, O., and Lee, Y. C. (1988) *Anal. Biochem.* **174**, 459–470.
37. Townsend, R. R., Hardy, M. R., Cumming, D. A., Carver, J. P., and Bendiak, B. (1989) *Anal. Biochem.* **182**, 1–8.
38. Hardy, M. R., and Townsend, R. R. (1988) *Proc. Natl. Acad. Sci. U.S.A.* **85**, 3289–3293.
39. Warren, L., and Glick, M. C. (1966) in *Methods in Enzymology* (Eastabrook, R. W., and Pullman, M. E., Eds.), Vol. 8, pp. 131–133, Academic Press, San Diego.
40. Schauer, R. (1978) in *Methods in Enzymology* (Ginsburg, V., Ed.), Vol. 50, pp. 64–89, Academic Press, San Diego.
41. Jourdian, G. W., Dean, L., and Roseman, S. (1971) *J. Biol. Chem.* **246**, 430–435.
42. Mellis, S. J., and Baenziger, J. U. (1981) *Anal. Biochem.* **114**, 276–280.
43. Bergh, M. L., Koppen, P., and Van den Eijnden, D. H. (1981) *Carbohydr. Res.* **94**, 225–229.
44. Mellis, S. J., and Baenziger, J. U. (1983) *Anal. Biochem.* **134**, 442–449.
45. Mellis, S. J., and Baenziger, J. U. (1983) *J. Biol. Chem.* **258**, 11,546–11,556.
46. Baenziger, J. U., and Natowicz, M. (1981) *Anal. Biochem.* **112**, 357–361.
47. Gross, H. J., Bunsch, A., Paulson, J. C., and Brossmer, R. (1987) *Eur. J. Biochem.* **168**, 595–602.
48. Higa, H. H., and Paulson, J. C. (1985) *J. Biol. Chem.* **260**, 8838–8849.
49. Kamerling, J. P., Schauer, R., Shukla, A. K., Stoll, S., van Halbeek, H., and Vliegthart, J. F. (1987) *Eur. J. Biochem.* **162**, 601–607.
50. Shukla, A. K., Schauer, R., Unger, F. M., Zahringer, U., Rietschel, E. T., and Brade, H. (1985) *Carbohydr. Res.* **140**, 1–8.
51. Shukla, A. K., Schauer, R., Schade, U., Moll, H., and Rietschel, E. T. (1985) *J. Chromatogr.* **337**, 231–238.
52. Herrler, G., Rott, R., Klenk, H. D., Muller, H. P., Shukla, A. K., and Schauer, R. (1985) *EMBO J.* **4**, 1503–1506.