

## Instructions for ATAC-seq sample submission for UCSD Investigators

Please contact <u>epigenome@ucsd.edu</u> to discuss your project with Center for Epigenomics staff before submitting samples. Once you have spoken to Center staff and are ready to submit samples, please:

(1) Ensure that your samples are prepared according to guidelines <u>here</u>.

(2) Complete sample manifest (download blank copy here).

- Please be sure to complete all fields in "Samples" tab.
- Please be sure to review policies and affirm consent on "Policies" tab.
- Please be sure that your samples are in 1.5 mL tubes clearly labeled with the same Sample ID as indicated on your sample submission form.
- Please also label side of tubes with your last name.
- (3) Send manifest to <u>epigenome@ucsd.edu</u> with the subject line: "ATAC-seq sample manifest - [your last name here].

(4) Please Drop samples off Monday-Thursday between 10am-4pm at our drop-off freezer (map here).